




## TECHNICAL DATA SHEET

Product Name:	Disposable Single Use Angiography Pack
Item No:	ANG-STP0001
Classification:	A
Description:	Majority components are made of SMMS, by cutting, punching and processing with certain protective characteristics, including 7 accessories: table cover, Mayo stand cover, reinforced surgical gown, C-section drape, baby drape, etc.
Intended Use:	Clinical surgery in relevant departments of medical institutions.
Material:	65gsm absorbent PP+PE+PP
Manufacturer / Assembler:	Graceheal Medic Sdn. Bhd
Country of Origin:	China / Malaysia
Specification:	





Component	Description	Material	Color	Quantity
Table Cover	140x190cm (60x190cm)	PP + PE	Blue	1
Mayo Stand Cover	76x145cm (48x60cm)	PP + PE	Blue	1
Surgical Gown	XL 130x150cm	SMMS	Blue	2
Angiography Drape	310x380cm (140x140cm) with absorbent pad and four circular fenestration	SMMS	Blue	1
Spong Brush	-	Sponge	White	1
Plastic Round Basin	Φ13cmx6cm	Plastic	Blue	1
Gauze	10 x 10cm	-	White	10
Hand Towel	30x40cm	Paper	White	2

**Packaging:**

- 1) Carton Requirements: Single layer
- 2) Packing Quantity: 1 set /header pouch, 6 set /ctn

**Storage:**

- 1) Store in dry, clean conditions in original packaging.
- 2) Avoid direct sunlight, high temperature or humidity.
- 3) Storage temperature ranging from -5°C to +45°C and with relative humidity below 80%.

**Shelf Life:**

36 months from the date of manufacture when stored as stated above.



**GRACEHEAL MEDIC SDN.BHD.**  
202301020220 ( 1514142-K )

23-3 A, JALAN PJU 1/3D,  
SUNWAY MAS COMMERCIAL CENTRE  
47301 PETALING JAYA,  
SELANGOR DARUL EHSAN,  
MALAYSIA.

+ 6014 333 8382 | + 6016 655 1425  
info@gracehealmedic.com

**Cautions and disposal:**

- 1) Do not contact with heavy oils, sparks or flames.
- 2) Do not use in the environment with conditions of excessive heat.
- 3) Do not use if damaged, heavily contaminated or in accordance with local practice.
- 4) Dispose the contaminated products with care and in accordance with national regulations.



Keep away from sunlight



Medical Device



Keep dry



Do not reuse



Temperature limit



Don't use when packing damaged  
and consult instructions for use



EO Sterile